



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Needle holder, vacuum tubes, box/100

**Product Code :**

LBNY-0001-21700016



**Description :**

Blood collection needle holder for use with vacuum tubes, and vacuum tube needles.  
For use in clinical phlebotomy.  
Material; clear polypropylene plastic  
Fits all standard, 13mm, and 16mm diameter vacuum tubes, and all common needle units.  
Disposable  
Blood collection needle holder.

**Technical Specification :**

Shelf life:  
36 months.

Weight and Volume:  
Estimated volume: 0.0073 m3.  
Estimated weight: 0.5 kg.



**Laboratory Instrument India**

