



Email : sales@laboratoryinstrumentindia.com

Product Name :

Swab, cotton-tip, sterile tube, box/100

Product Code :

LBNY-0001-21700017



Description :

Used for various clinical and laboratory purposes.

Tube material: Polypropylene.

Plastic handle stick.

Sterile, individually packed in a tube.

Rayon/cotton tipped.

Length:12cm.

Technical Specification :

Shelf life:

36 months.

Weight and Volume:

Estimated volume: 0.00035 m3.

Estimated weight: 0.058 kg.



Laboratory Instrument India