



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Container, sample, 50 ml, box/100

**Product Code :**  
LBNY-0001-21700018



**Description :**

Capacity Approx. 50ml.  
Non-sterile, Single use/disposable.  
Screw cap, wide mouth, leakproof.  
Clear plastic.  
Sample container used to collect a wide range of clinical specimens.

**Technical Specification :**

Shelf life:  
60 months

Estimated volume: 0.014 m3.



**Laboratory Instrument India**