



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Box,infectious sample transport

**Product Code :**  
LBNY-0001-21700019



**Description :**

Transportation of Category-A infectious substance, with approved triple packaging.

Leakproof box, with inserts providing barriers.

Compliant design for transportation (infectious substance affecting humans, category A), and approved for Class 6.2 hazardous materials.

Packing capacity:

Tertiary: Very strong cardboard (sturdy enough to withstand a fall from 3 meters high), with minimal dimensions (to reduce air-freight cost), includes integrated space for 0.3L cool-pack (for short-distance transportation) or dry-ice (for long-distance transportation).

Primary pack: Loads up to six standard 5-7ml vacuum specimen tubes.

Secondary pack: Plastic container, Approx. 60x160mm (diam x h), fitted with absorbent material.

**Technical Specification :**

Weight and Volume:

Estimated volume: 0.014 m3.

Estimated weight: 2 kg.



**Laboratory Instrument India**

