



Email : sales@laboratoryinstrumentindia.com

Product Name :

Tube,blood coll,w/capil,0.2 ml, EDTA 100

Product Code :

LBNY-0001-21700020



Description :

Used for capillary blood collection for clinical Lab tests.

Capacity: 0.2 ml.

Material: plastic.

Capped with a pierceable flat vacuum seal.

Indicates maximum fill with clearly visible line.

Blood collection tube

Prefilled with K2 EDTA,

Capillary tube included.

Technical Specification :

Shelf life:

18 Months.

Weight and Volume:

Estimated volume: 0.0027m³.

Estimated weight: 0.32 kg.

Laboratory Instrument India

