




Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube,blood coll,w/capil,0.2 ml, EDTA 100

**Product Code :**

LBNY-0001-21700020



**Description :**

Used for capillary blood collection for clinical Lab tests.

Capacity: 0.2 ml.

Material: plastic.

Capped with a pierceable flat vacuum seal.

Indicates maximum fill with clearly visible line.

Blood collection tube

Prefilled with K2 EDTA,

Capillary tube included.

**Technical Specification :**

Shelf life:

18 Months.

Weight and Volume:

Estimated volume: 0.0027m<sup>3</sup>.

Estimated weight: 0.32 kg.

---

# Laboratory Instrument India

