



Email : sales@laboratoryinstrumentindia.com

Product Name :

Container, sharps, leak-proof plastic,4L

Product Code :

LBNY-0001-21700021



Description :

Made of rigid plastic.
Capacity up to 4.7 L.
Clear lid to view fill level.
Leak proof, puncture resistant.
Built-in needle removal port.
Upright and stable during use.
Properly labelled with a biohazard label.
Able to close with a tight-fitting, puncture-resistant lid without sharps being able to come out.
Leak-proof, puncture-resistant plastic container for disposing potentially infectious sharps.

Technical Specification :

Weight and Volume:
Estimated volume: 0.0059 m3.
Estimated weight: 0.4 kg.



Laboratory Instrument India

