




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Rack, drying DBS cards, 10 pos, box/10

**Product Code :**  
LBNY-0001-2170022



**Description :**

Material: Cardboard.  
Sturdy base provides optimal stability.  
Self-supporting vertical rack for drying DBS cards.  
Provides horizontal position for 10 cards.  
Used to air dry the blood spot samples collected for early infant diagnosis of HIV.  
Single use.

**Technical Specification :**

Weight and Volume:  
Estimated volume:0.0059 m3.  
Estimated weight: 0.4 kg.



**Laboratory Instrument India**