




Email : sales@laboratoryinstrumentindia.com

Product Name :
Rack, drying DBS cards, 10 pos, box/10

Product Code :
LBNY-0001-2170022



Description :

Material: Cardboard.
Sturdy base provides optimal stability.
Self-supporting vertical rack for drying DBS cards.
Provides horizontal position for 10 cards.
Used to air dry the blood spot samples collected for early infant diagnosis of HIV.
Single use.

Technical Specification :

Weight and Volume:
Estimated volume:0.0059 m3.
Estimated weight: 0.4 kg.



Laboratory Instrument India