




Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Rack, drying DBS cards, 10 pos, box/10

**Product Code :**

LBNY-0001-2170022



**Description :**

Material: Cardboard.

Sturdy base provides optimal stability.

Self-supporting vertical rack for drying DBS cards.

Provides horizontal position for 10 cards.

Used to air dry the blood spot samples collected for early infant diagnosis of HIV.

Single use.

**Technical Specification :**

Weight and Volume:

Estimated volume:0.0059 m3.

Estimated weight: 0.4 kg.



**Laboratory Instrument India**