




Email : sales@laboratoryinstrumentindia.com

Product Name :

Lancet, safety, 2.0 mm, disp, BOX/100

Product Code :

LBNY-0001-21700023



Description :

Used for capillary blood collection
In-built mechanism to retract needle tip preventing reuse or injury.
Instant trigger launches spring-based blade.
Sterile, disposable.
Blade depth approx. 2.0mm.

Technical Specification :

Shelf life:
60 months.

Weight and Volume:
Estimated volume:0.0015 m3.
Estimated weight:0.3 kg.



Laboratory Instrument India

