




Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Lancet, safety, 2.0 mm, disp, BOX/100

**Product Code :**

LBNY-0001-21700023



**Description :**

Used for capillary blood collection  
In-built mechanism to retract needle tip preventing reuse or injury.  
Instant trigger launches spring-based blade.  
Sterile, disposable.  
Blade depth approx. 2.0mm.

**Technical Specification :**

Shelf life:  
60 months.

Weight and Volume:  
Estimated volume:0.0015 m3.  
Estimated weight:0.3 kg.



**Laboratory Instrument India**

