



Email : sales@laboratoryinstrumentindia.com

Product Name :

Lancet, safety, 2.4 mm, disp, BOX/100

Product Code :

LBNY-0001-21700024



Description :

Instant trigger launches spring-based blade
In-built mechanism to retract needle tip preventing reuse or injury.
Used for capillary blood collection.
Blade depth approx. 2.4mm
Single-use disposable
Sterile.

Technical Specification :

Shelf life:
60 months.

Weight and Volume:
Estimated volume: 0.0015 m3.
Estimated weight: 0.3 kg.



Laboratory Instrument India

