



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

Lancet, safety, 2.4 mm, disp, BOX/100

**Product Code :**

LBNY-0001-21700024



**Description :**

Instant trigger launches spring-based blade  
In-built mechanism to retract needle tip preventing reuse or injury.  
Used for capillary blood collection.  
Blade depth approx. 2.4mm  
Single-use disposable  
Sterile.

**Technical Specification :**

Shelf life:  
60 months.

Weight and Volume:  
Estimated volume: 0.0015 m3.  
Estimated weight: 0.3 kg.



**Laboratory Instrument India**

