



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Box,infectious sample transport

**Product Code :**  
LBNY-0001-21700025



**Description :**

Primary receptacle: Leakproof box, with inserts providing barriers.  
Secondary packing: plastic container, can hold 4 tubes up to 20ml. Approx. 40x150mm (diam x h), fitted with absorbent material.  
Transport of Category-B infectious substance, with approved triple packaging  
Compliant design for transportation.

**Technical Specification :**

Weight and Volume:  
Estimated volume: 0.014 m3.  
Estimated weight: 2 kg.



**Laboratory Instrument India**