



Email : sales@laboratoryinstrumentindia.com

Product Name :
Box,infectious sample transport

Product Code :
LBNY-0001-21700025



Description :

Primary receptacle: Leakproof box, with inserts providing barriers.
Secondary packing: plastic container, can hold 4 tubes up to 20ml. Approx. 40x150mm (diam x h), fitted with absorbent material.
Transport of Category-B infectious substance, with approved triple packaging
Compliant design for transportation.

Technical Specification :

Weight and Volume:
Estimated volume: 0.014 m3.
Estimated weight: 2 kg.



Laboratory Instrument India