



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Tube, vacuum, EDTA, gel,PPT,6ml, BOX-100

**Product Code :**

LBNY-0001-21700027



**Description :**

Plasma collection for clinical laboratory tests. Intended for the purposes of molecular diagnostic testing.

Material: clear nonfragile PET plastic and leak-proof rubber cap

Plastic shielded rubber cap to protect personnel from direct contact with blood.

Blood collection vacuum tube.

Volume 6m,

Sterile

Standard color-coded cap.

Pen/Marker writable label space.

Capped, with needle-pierceable vacuum seal cap.

Additive; Inert plasma separation gel & K2 EDTA anticoagulant.

**Technical Specification :**

Shelf life: 16 months.

Weight & volume:

Estimated volume: 0.00495 m3.

Estimated weight: 0.91 kg.

---

# Laboratory Instrument India

