



Email : sales@laboratoryinstrumentindia.com

Product Name :
Blood grouping,anti-A,IgM,10ml,box/10

Product Code :
LBNY-0001-2180001



Description :

Anti-A blood grouping reagent
Volume: vial of 10ml
Dropper bottle
Contains monoclonal murine IgM antibody
Box of 10 vials.

Instructions for use:
Used for blood transfusion services
Storage and transportation: between 2 and 8C.
Used in wide range of laboratory services from basic laboratory to reference centre.
Shelf life: 36 months after manufacturing date.

Technical Specification :

Weight - Volume
Estimated volume: 0.0008m3
Estimated weight: 0.34kg

Laboratory Instrument India

