



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
SMD Rework Station (Hot Air)

**Product Code :**  
LBNY-0002-610016



**Description :**

SMD Rework Station (Hot Air)

**Technical Specification :**

MCU and PID algorithm controlled temperature be more accurate.  
Soldering iron and hot air station integrated to save space.  
Dual High definition LCD display to show various parameters make operation convenient.  
Digital temperature calibration and switchable temperature units make it versatile.  
Diaphragm pump with high pressure make it suitable for various nozzles.  
Intelligence cooling system and deferred power-off function to extend hand piece lifespan.  
Patented heater for obtaining fast heating-up and extending lifespan.  
Three preset buttons to switch between different temperature and air grades combinations.  
Automatic standby function to save more energy.  
Intelligent malfunction prompt to let user locate causes fast and accurately.  
Control warm (Temperature Adjustable)  
Handpiece Length: 196MM . 110V 220V available anti-static.

SMD Rework Station (Hot Air), SMD Rework Station (Hot Air) Equipments, SMD Rework Station (Hot Air) Tools, SMD Rework Station (Hot Air) Tool Kits, SMD Rework Station (Hot Air) Manufacturers, SMD Rework Station (Hot Air) Suppliers from India, China, Kenya



**Laboratory Instrument India**