



Email : sales@laboratoryinstrumentindia.com

Product Name :
Alcohol Thermometer

Product Code :
LBNY-0003-830001



Description :

Alcohol Thermometer, Alcohol Thermometer Equipments, Alcohol Thermometer Tools, Alcohol Thermometer Tool Kits, Alcohol Thermometer Manufacturers, Alcohol Thermometer Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

Alcohol Thermometer
Measurement range: 0 - 150°C
Scale: 1°C
Spirit filled (red).



Laboratory Instrument India

