



Email : sales@laboratoryinstrumentindia.com

Product Name :

T-shirt, Child 4-6yrs,cols/PAC-4

Product Code :

LBNY-0001-21900011



Description :

Children's T-shirt, assorted colours, 4-6 years, pack of 4.

Technical Specification :

Material Composition: 45% cotton 45% polyester and 10% elastin. 165 gsm, flat knit with a round neck and short sleeves. The dimensions of the neck allow comfortable wearing and removal by children.

Children's t-shirt, short sleeves. Round neck and twin needle hems, with a hidden seam at the neck.

The material is colour-fast, durable and easy to wash. The material and design are suitable for extensive wear in strong sunlight – the colour is durable and does not bleach easily in strong sunlight. The colour (material, thread, and print) is able to withstand regular and intensive washing cycles with a mild washing powder. There is no leaching of colour when the t-shirt comes into contact with water.



Laboratory Instrument India