



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

T-shirt, Child 8-10yrs,cols/PAC-4

**Product Code :**

LBNY-0001-21900012



**Description :**

Material Composition: 45% cotton 45% polyester and 10% elastin. 165 gsm, flat knit with a round neck and short sleeves. The dimensions of the neck allow comfortable wearing and removal by children.

Children's t-shirt, short sleeves. Round neck and twin needle hems, with a hidden seam at the neck.

The material is colour-fast, durable and easy to wash. There is no leaching of colour when the t-shirt comes into contact with water. The colour (material, thread, and print) is able to withstand regular and intensive washing cycles with a mild washing powder. The material and design are suitable for extensive wear in strong sunlight – the colour is durable and does not bleach easily in strong sunlight.

**Technical Specification :**

T-shirt, Child 8-10yrs,cols/PAC-4



**Laboratory Instrument India**