



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**

T-shirt, Child 12-14yrs,cols/PAC-4

**Product Code :**

LBNY-0001-21900013



**Description :**

Children's T-shirt, assorted colours, 12-14years, pack of 4.

**Technical Specification :**

Material Composition: 45% cotton 45% polyester and 10% elastin. 165 gsm, flat knit with a round neck and short sleeves. The dimensions of the neck allow comfortable wearing and removal by children.

Children's t-shirt, short sleeves. Round neck and twin needle hems, with a hidden seam at the neck.

The material is colour-fast, durable and easy to wash. The colour (material, thread, and print) is able to withstand regular and intensive washing cycles with a mild washing powder. The material and design are suitable for extensive wear in strong sunlight – the colour is durable and does not bleach easily in strong sunlight. There is no leaching of colour when the t-shirt comes into contact with water.



**Laboratory Instrument India**