



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Aquarium with Filter and Air Pump

**Product Code :**  
LBNY-0003-830002



**Description :**

Aquarium with Filter and Air Pump

**Technical Specification :**

Aquarium with filter and air pump

There are four suction cups on the back to hold the filter on the tank wall Made of high quality sturdy ABS housing, stainless steel shaft, durable, fresh water and sea water. Features such as filtration, oxygen absorption, pumping, etc.

Removable filter cartridge for better filtration. Small size, easy to hide Duckbill nozzle design, support adjustment of the angle of the water outlet and airflow direction, can be adjusted by 180 Multi-layer biochemical filtration system and large-capacity filtration space for enhanced filtration.

Capacity: 50L

Size: 50 x 25x 30 cm

Water filter filtering

Capacity: 100 l/h

Air pump-110-130 l/h



**Laboratory Instrument India**

