



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cold box, B.Kings

**Product Code :**  
LBNY-0001-210001



**Description :**

Vaccine storage capacity: 7 litres  
Cold box, Small, Short range  
Cold life at +43°C: 89 hours

**Technical Specification :**

Vaccine storage capacity: 7 litres  
Vaccine cold box, Small, Short range  
Cold life at 43°C without openings: 89 hours  
With hinged lid and carrying handles.  
Weight fully loaded: 23 Kg  
Supplied with 24 x 0.4L icepacks as a standard.

**Dimensions:**

Internal dimensions (LxWxH cm): 35 x 28.5 x 27  
External dimensions (LxWxH cm): 49 x 44 x 39.5



**Laboratory Instrument India**

