



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box, B.Kings

Product Code :
LBNY-0001-210001



Description :

Vaccine storage capacity: 7 litres
Cold box, Small, Short range
Cold life at +43°C: 89 hours

Technical Specification :

Vaccine storage capacity: 7 litres
Vaccine cold box, Small, Short range
Cold life at 43°C without openings: 89 hours
With hinged lid and carrying handles.
Weight fully loaded: 23 Kg
Supplied with 24 x 0.4L icepacks as a standard.

Dimensions:

Internal dimensions (LxWxH cm): 35 x 28.5 x 27
External dimensions (LxWxH cm): 49 x 44 x 39.5



Laboratory Instrument India

