




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Vacc.carr.,B.Kings

**Product Code :**  
LBNY-0001-210002



**Description :**

Vaccine storage capacity: 1.7 litres  
Vaccine carrier, Large, Short range  
Cold life at +43°C without openings: 38 hours  
Cold life +43°C: 38 hours  
With lid and shoulder strap.  
Supplied with 4 x 0.4L icepacks as a standard.  
Weight fully loaded: 4 Kg.

**Technical Specification :**

Dimensions :  
Internal dimensions (LxWxH cm): 17 x 17 x 18.5  
External dimensions (LxWxH cm): 26 x 25 x 28.5



**Laboratory Instrument India**