




Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc.carr.,B.Kings

Product Code :
LBNY-0001-210002



Description :

Vaccine storage capacity: 1.7 litres
Vaccine carrier, Large, Short range
Cold life at +43°C without openings: 38 hours
Cold life +43°C: 38 hours
With lid and shoulder strap.
Supplied with 4 x 0.4L icepacks as a standard.
Weight fully loaded: 4 Kg.

Technical Specification :

Dimensions :
Internal dimensions (LxWxH cm): 17 x 17 x 18.5
External dimensions (LxWxH cm): 26 x 25 x 28.5



Laboratory Instrument India