



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vacc Carrier Leff

**Product Code :**  
LBNY-0001-210003



**Description :**

Vaccine storage capacity: 1.7 liters  
With lid, shoulder strap and backpack.  
Cold life 33 hours 42 minutes at 43°C without openings.  
Supplied with 4 x 0.6 liter icepacks WP0.6L as a standard.  
Weight fully loaded: 8 kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH): 10 x 10 x 19 cm  
External dimensions (LxWxH): 30.8 x 30.8 x 30 cm



**Laboratory Instrument India**