



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc Carrier Leff

Product Code :
LBNY-0001-210003



Description :

Vaccine storage capacity: 1.7 liters
With lid, shoulder strap and backpack.
Cold life 33 hours 42 minutes at 43°C without openings.
Supplied with 4 x 0.6 liter icepacks WP0.6L as a standard.
Weight fully loaded: 8 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH): 10 x 10 x 19 cm
External dimensions (LxWxH): 30.8 x 30.8 x 30 cm



Laboratory Instrument India