



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vac.carr.,B.Kings

Product Code :
LBNY-0001-210004



Description :

Vaccine storage capacity: 2.6 litres
Vaccine carrier, Extra-large, Long range
With lid and shoulder strap.
Supplied with 4 x 0.6L icepacks as a standard.
Cold life at +43°C without openings: 43 hours.
Weight fully loaded 4.5 kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 19 x 19 x 22.5
External dimensions (LxWxH cm): 26 x 26 x 32



Laboratory Instrument India