



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carr,BBVC-23

Product Code :
LBNY-0001-210005



Description :

Vaccine storage capacity: 0.9 litres
Vaccine carrier, Small, Short range
Cold life at 4°C without openings: 17 hours 50 minutes
With lid and shoulder strap.
Supplied with 2 x 0.3L icepacks as a standard.
Weight fully loaded: 2.4 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 17.4 x 10.5 x 12.0
External dimensions (LxWxH cm): 24.6 x 18.0 x 21.5



Laboratory Instrument India