



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carrier

Product Code :
LBNY-0001-210006



Description :

Vaccine storage capacity: 2.5 liters
Vaccine carrier, Extra Large, Long range
With lid and shoulder strap.
Cold life at 43°C without openings: 50 hours and 12 minutes
Supplied with 4 x 0.6L icepacks as a standard.
Weight fully loaded 6.36 kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 12.20 x 12.20 x 19.00
External dimensions (LxWxH cm): 27.00 x 27.00 x 32.00



Laboratory Instrument India