



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold Box AICB

Product Code :
LBNY-0001-210007



Description :

Vaccine storage capacity: 8 litres
Vacc. carrier Cold box, Small, Short range
With hinged lid and carrying handles.
Cold life at 43°C without openings: 84 hours 7 minutes
Supplied with 15 x 0.3L icepacks as a standard.
Weight fully loaded: 22.41 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 34.6 x 28 x 24
External dimensions (LxWxH cm): 54.4 x 44.5 x 42



Laboratory Instrument India