



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold Box Long term storage 5.4L

Product Code :
LBNY-0001-210008



Description :

Cold box, Long term vaccine storage.
Storage capacity: 5.4L
Operating Ambient Temperature: +43°C
Vaccine storage Temperature: +2°C to + 8°C
Cold life at +43°C: 35 Days.

Technical Specification :

Ice packs should be placed in a deep freezer, frozen between -5°C to -20°C for at least 48 hours and then conditioned as explained in the User Manual. Then these icepacks can be moved into the Cold box to cool the device to +2°C to +8°C.

Cold life: Cold box can keep vaccines between +2°C to +8°C temperature for 35-days if properly used at 43°C.
With storage capacity: 5.4L
With Vacuum Dewar insulation
Internal: Diam 31.1 x 33.2 cm
External: Diam 52.8 x 74.7 cm
Coolant pack size: 274*101*74 mm
Coolant pack Type: 1L ice packs
Empty/Full weight: 22kg/39.5kg
Number of water packs used with the cold box: 8 pcs of ice packs (ice blocks)/every equipment.



Laboratory Instrument India