



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cold box 20 litres

**Product Code :**  
LBNY-0001-210009



**Description :**

Vaccine storage capacity: 20 litres  
Cold box, Large, Long range  
With hinged lid and carrying handles.  
Cold life at 43°C without openings: 134.6 hours  
Supplied with 24 x 0.6L icepacks, as a standard.  
Weight fully loaded: 38.9 Kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 49.6 x 33.4 x 26.4  
External dimensions (LxWxH cm): 71 x 55 x 49.9



**Laboratory Instrument India**