



Email : sales@laboratoryinstrumentindia.com

Product Name :
Binocular Microscope

Product Code :
LBNY-0003-830003



Description :

Binocular Microscope

Technical Specification :

Binocular Microscope

Brightfield, binocular with double layer stage; objectives 4x, 10x, 40x, 100x

Illumination: Light source type: white LED; with light intensity control

Focusing: Coaxial coarse and fine focusing mechanism (graduated, 0.002 mm) with upper stop, to prevent the contact between objective and specimen. Adjustable tension of coarse focusing knob.

Stage: Double layer with mechanical sliding stage, size 125 x 116 mm, X-Y movement range 76 x 30. Specimen holder for one slide. Vernier scale on the two axes, accuracy 0,1 mm.

Nosepiece: Quadruple revolving nosepiece

Head Binocular observation head, inclined 30° and rotatable 360°. Diopter adjustment on left eyepiece.

Interpupillary adjustment 48 - 75 mm

Eyepieces. Widefield eyepieces WF10X/18 with field number 18 mm.

Condenser: Abbe condenser, N.A. 1.25 precentered.

Instruction manual and dust cover is included.

Include:

The Eyepiece reticule- 10x

Stage micrometer.



Laboratory Instrument India