



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box 7 litres

Product Code :
LBNY-0001-2100010



Description :

Vaccine storage capacity: 7 litres
Cold box, Small, Long range
Cold life at 43°C without openings: 114.9 hours
With hinged lid and carrying handles.
Supplied with 14 x 0.6L icepacks as a standard.
Weight fully loaded: 23.3 Kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 34 x 26 x 27
External dimensions (LxWxH cm): 55 x 47.5 x 49.9



Laboratory Instrument India