



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Cold box, Nilkamal

**Product Code :**  
LBNY-0001-2100011



**Description :**

Vaccine storage capacity: 16 litres  
Cold box, Large, Short range  
Cold life at 43°C without openings: 73 hours 32 minutes  
With hinged lid and carrying handles.  
Supplied with 24 x 0.6L icepacks as a standard.  
Weight fully loaded: 34.8 Kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 48 x 48 x 21  
External dimensions (LxWxH cm): 65 x 65 x 37



**Laboratory Instrument India**