



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cold box 18 litres

**Product Code :**  
LBNY-0001-2100012



**Description :**

Vaccine storage capacity: 18 litres  
Cold box, Large, Long range  
With hinged lid and carrying handles.  
Cold life at 43°C without openings: 140 hours  
Supplied with 44 x 0.4L icepacks as a standard.  
Weight fully loaded: 48.87 Kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 52 x 37 x 28.2  
External dimensions (LxWxH cm): 76.1 x 61.1 x 51.3



**Laboratory Instrument India**