



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cold box, AOV

**Product Code :**  
LBNY-0001-2100013



**Description :**

Vaccine storage capacity: 6 litres  
Cold box, Small, Long range  
With hinged lid and carrying handles.  
Cold life +43°C without opening: 132 hrs 18 mins  
Supplied with 26 x 0.4 litre icepacks as a standard.  
Weight fully loaded: 25.40 Kg.

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 40.30 x 28.20 x 21.10  
External dimensions (LxWxH cm): 62.40 x 50.20 x 42.60



**Laboratory Instrument India**