



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box, AOV

Product Code :
LBNY-0001-2100013



Description :

Vaccine storage capacity: 6 litres
Cold box, Small, Long range
With hinged lid and carrying handles.
Cold life +43°C without opening: 132 hrs 18 mins
Supplied with 26 x 0.4 litre icepacks as a standard.
Weight fully loaded: 25.40 Kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 40.30 x 28.20 x 21.10
External dimensions (LxWxH cm): 62.40 x 50.20 x 42.60



Laboratory Instrument India