



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carrier,AOV

Product Code :
LBNY-0001-2100014



Description :

Vaccine storage capacity: 0.8 litres
Vaccine carrier, Small, Short range
With lid and shoulder strap.
Cold life at 43°C without openings: 20 hours 16 minutes
Supplied with 2 x 0.4L icepacks as a standard.
Weight fully loaded: 2.18 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 17.30 x 10.30 x 12.10
External dimensions (LxWxH cm): 25.00 x 18.00 x 21.00



Laboratory Instrument India