



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cold box 5.5 litres

**Product Code :**  
LBNY-0001-2100015



**Description :**

Cold box, Small, Long range  
Vaccine storage capacity: 5.5 litres  
With hinged lid and carrying handles.  
Cold life at 43°C without openings: 107.5 hours  
Supplied with 15 x 0.6L icepacks as a standard.  
Weight fully loaded: 22.33 Kg

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 35 x 28 x 24  
External dimensions (LxWxH cm): 54 x 44.5 x 41.5



**Laboratory Instrument India**