



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box 5.5 litres

Product Code :
LBNY-0001-2100015



Description :

Cold box, Small, Long range
Vaccine storage capacity: 5.5 litres
With hinged lid and carrying handles.
Cold life at 43°C without openings: 107.5 hours
Supplied with 15 x 0.6L icepacks as a standard.
Weight fully loaded: 22.33 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 35 x 28 x 24
External dimensions (LxWxH cm): 54 x 44.5 x 41.5



Laboratory Instrument India