



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cold box 22.5 litres

Product Code :
LBNY-0001-2100016



Description :

Vaccine storage capacity 22.5 litres
Cold box, Large, Long range
With hinged lid and carrying handles.
Cold life at 43°C without openings: 128 hours 16 minutes
Supplied with 50 x 0.3L icepacks as a standard.
Weight fully loaded: 48.02 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 51.8 x 37 x 27.7
External dimensions (LxWxH cm): 76.5 x 61.2 x 51.5



Laboratory Instrument India