



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carrier 2.9 litres

Product Code :
LBNY-0001-2100017



Description :

Vaccine storage capacity: 2.9 litres
Vaccine carrier, Extra-large, Long range
With lid and shoulder strap.
Cold life at 43°C without openings: 51 hours 8 minutes
Supplied with 4 x 0.6L icepacks as a standard.
Weight fully loaded: 5.36 Kg.

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 20 x 20 x 19
External dimensions (LxWxH cm): 29 x 29 x 32.7



Laboratory Instrument India