



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vacc. carrier,BMed

Product Code :
LBNY-0001-2100018



Description :

Vaccine storage capacity: 3 litres
Vaccine carrier, Extra-large, Long range
With lid and shoulder strap.
Cold life at 43°C without openings: 30.3 hours
Supplied with 6 x 0.3L & 1 x 0.6L icepacks as a standard.
Weight fully loaded: 7.3 Kg

Technical Specification :

Dimensions:
Internal dimensions (LxWxH cm): 26.0 x 15.6 x 18.6
External dimensions (LxWxH cm): 36.2 x 28.3 x 29.9



Laboratory Instrument India