



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Vacc. carrier,BMed

**Product Code :**  
LBNY-0001-2100018



**Description :**

Vaccine storage capacity: 3 litres  
Vaccine carrier, Extra-large, Long range  
With lid and shoulder strap.  
Cold life at 43°C without openings: 30.3 hours  
Supplied with 6 x 0.3L & 1 x 0.6L icepacks as a standard.  
Weight fully loaded: 7.3 Kg

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 26.0 x 15.6 x 18.6  
External dimensions (LxWxH cm): 36.2 x 28.3 x 29.9



**Laboratory Instrument India**