



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cold box, BMed

**Product Code :**  
LBNY-0001-2100019



**Description :**

Cold box, Small, Short range  
Vaccine storage capacity: 6 litres  
With hinged lid and carrying handles.  
Cold life at 43°C without openings: 57.9 hours  
Supplied with 2 x 0.3L and 10 x 0.6L icepacks as a standard.  
Weight fully loaded: 16.4 Kg

**Technical Specification :**

Dimensions:  
Internal dimensions (LxWxH cm): 46 x 18 x 24.5  
External dimensions (LxWxH cm): 58.8 x 28.8 x 43.7



**Laboratory Instrument India**